# VANTAstar™

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


Injectors enable reagent delivery to any plate format from 6 - 384 wells.

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- Standard Curve Wizard for a step-by-step standard curve calculation
- Automatic DNA/RNA concentration determination
- Data display as bar charts, box plots, violin plots etc.
- Spectral view and analysis
- Background and baseline correction
- Signal interpolation: linear or cubic spline
- Various curve fit models including linear, 4-/5-parameter, polynomial and user-defined fit models
- Enzyme kinetic analysis using various models
- EC<sub>50</sub> calculation with confidence intervals
- Binding rates and constants determination
- ANOVA, Student's t-test or multiple comparisons
- Performance evaluation: signal-to-blank, signal-to-noise, %CV, Z-prime, etc.
- Automatic data processing using predefined templates

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
- DNA, RNA, and protein quantification
- ELISA and DELFIA immuno-assays
- Bacterial growth (OD<sub>600</sub>)
- Cell viability/toxicity
- Real-time cell-based assays
- Enzymatic activity
- Reporter gene assays
- And much more ...

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LANCE and DELFIA are registered trademarks of PerkinElmer, Inc.  
HTRF is a registered trademark of Cisbio Bioassays.  
LanthaScreen is a registered trademarks of Thermo Fisher Scientific.  
Transcreener is a registered trademark of Bellbrook Labs.  
Dual Luciferase Reporter is a trademark of Promega Corp.  
Mycoalert is a trademark of Lonza.  
THUNDER is a trademark of Bioauxilium.

VANTASTAR can include all the following functions or any combination of functions, and can also be upgraded with additional functions at any time in the future.

Detection modes	UV/vis absorbance Fluorescence intensity (incl. FRET) Luminescence (flash and glow) - incl. BRET Time-resolved fluorescence TR-FRET Fluorescence polarization	
Measurement modes	Top and bottom reading Endpoint and kinetic Sequential multi-excitation Sequential multi-emission Spectral scanning (absorbance, fluorescence, luminescence) Ratiometric measurements Well scanning	
Microplate formats	6- to 384-well plates, user-definable LVis Plate with 16 low volume microspots (2 µL)	
Microplate carrier	Robot compatible	
Light source	High energy xenon flash lamp	
Detectors	Low-noise photomultiplier tube CCD spectrometer	
Wavelength selection	Dual Linear Variable Filter (LVF) Monochromators™ Linear Variable Dichroic Mirror: separates excitation and emission LVF Monochromators Optical filters: excitation and emission slides hold up to 5 filters each LVF Monochromators + optical filters: use one for excitation and the other for emission UV/vis absorbance spectrometer: full spectra or 8 discrete wavelengths in <1 sec/well	
Optical filters	Excitation and emission slides for up to 4 filters each for monochromator units (VANTASTAR), and for up to 5 filters each for filter-only units (VANTASTAR F)	
Optical path	Top and bottom: free-air optical light path guided by motor-driven mirrors and dichroics	
Z-adjustment	Automatic focal height adjustment (0.1 mm resolution)	
Spectral range	<b>Filters</b>	FI, LUM, TRF, FP, TR-FRET: 240 - 740 nm
	<b>LVF Monochromators™</b>	FI, LUM: 320 - 740 nm
	<b>Linear Variable Dichroic</b>	340 - 740 nm
	<b>Spectrometer</b>	ABS: 220 - 1000 nm; wavelength precision: < 0.5 nm
Sensitivity	<b>FI filters (top)</b>	< 0.5 pM (< 10 amol/well FITC, 384sv, 20 µL)
	<b>FI filters (bottom)</b>	< 2.5 pM (< 125 amol/well FITC, 384g, 50 µL)
	<b>FI monochromator (top)</b>	< 0.8 pM (< 16 amol/well FITC, 384sv, 20 µL)
	<b>FI monochromator (bottom)</b>	< 3.5 pM (< 175 amol/well FITC, 384g, 50 µL)
	<b>FI dynamic range</b>	8 decades in a single measurement
	<b>LUM</b>	< 0.8 pM (< 15 amol/well ATP, 384sv white, 20 µL)
	<b>LUM dynamic range</b>	8 decades in a single measurement
	<b>TRF</b>	< 30 fM europium (< 2.4 amol/well, 384, 80 µL)
	<b>HTRF®</b>	Reader Control Kit [Eu] after 18h (384sv, 20 µL) Delta F > 700% (High Calibrator) Delta F > 25% (Low Calibrator)
	<b>FP</b>	< 1.2 mP SD at 1 nM FITC (384sv, 20 µL)
Read times	<b>ABS with spectrometer</b>	Selectable spectral resolution: 1, 2, 5, and 10 nm OD range: 0 to 4 OD; photometric resolution: 0.001 OD Accuracy: < 1% at 2 OD Precision: < 0.5% at 1 OD and < 0.8% at 2 OD Linearity: < 0.8% at 2.0 OD
	<b>Flying mode (1 flash)</b>	9 sec (96), 16 sec (384)
Shaking	Linear, circular, and double-circular with user-definable time and speed	
Incubation	+4°C above ambient up to 45°C or 60°C Minimum condensation concept: the upper heating plate of the incubation chamber operates at +0.5°C than the lower plate, heated optic system.	
Software	Multi-user Smart Control and MARS data analysis software included FDA 21 CFR Part 11 compliant Integrated fluorophore library	
Dimensions	Width: 35 cm, depth: 52 cm, height: 36 cm; weight: 27 kg	
<b>Optional accessories</b>		
Atmospheric Control Unit (ACU)	Actively regulates O <sub>2</sub> and CO <sub>2</sub> : 1-20%	
Reagent dispenser module	Up to 2 reagent injectors Individual injection volumes for each well: 3 to 500 µL (optionally up to 1 mL) Variable injection speed up to 420 µL/s Reagent back flushing Heater and stirrer plate Dimensions: width: 13 cm, depth: 31 cm, height: 26 cm; weight: 6 kg	
LVis Plate	Microplate designed to measure 16 low volume (2 µL) samples and standard cuvettes. Incorporating NIST-traceable filters and holmium oxide standards for instrument performance test. Sensitivity: < 2 ng/µL dsDNA	
THERMOstar	Microplate incubator and shaker	
Filters	Optimised for dyes, fluorophores and specific assays Filters for all applications from UV to NIR Customised filters available upon request	
Upgrades	Upgrades to include options such as additional detection modes, reagent injectors, extended temperature control, etc. are available. Please contact your local representative for more information.	

**Headquarters Germany**  
 BMG LABTECH GmbH  
 Allmendgrün 8  
 77799 Ortenberg  
 Tel. +49 781 96968 -0  
 sales@bmglabtech.com  
[www.bmglabtech.com](http://www.bmglabtech.com)